



# Recombinant *Oryza sativa* subsp. indica Homeobox-leucine zipper protein HOX2 (HOX2)

<b>Product Code</b>	CSB-YP771164OFF
<b>Abbreviation</b>	HOX2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q84U86
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. indica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMDLGLSLGL GLASQGSPTS STTTTSSPGA GSSSPWAAAL NSIVGDVRRD QAAAHAAAV GVGVGGEEMY QGRASTSPDS AAALSSASGK RERELERSGS GVDDDDGADG AGGRKKLRLS KDQAAVLEEC FKTHSTLNPK KQVALANRLG LRPRQVEVWF QNRRARTKLIK QTEVDCEYLK RWCERLADEN KRLEKELADL RALKAAPSPA SASAMQPSSS AAATLTMCPSS CRRVATAGAP HQPNHQQCHP KSNTTISSSS TAAA AVAVAG GNVLP SHCQF FPAAAAAADR TSQSTWNAAA PLVTRELF
<b>Source</b>	Yeast
<b>Target Names</b>	HOX2
<b>Protein Names</b>	Recommended name: Homeobox-leucine zipper protein HOX2 Alternative name(s): HD-ZIP protein HOX2 Homeodomain transcription factor HOX2 OsHox2
<b>Expression Region</b>	1-308
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.