



# Recombinant Rice dwarf virus Suppressor of RNA-mediated gene silencing

<b>Product Code</b>	CSB-YP771188REN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q85447
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rice dwarf virus (isolate Akita) (RDV)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEVDTATFVR LHELLCAHE GPSIISKFDA IKKVKLGTLA NQSGGVNNIT EAFLAKLRNF ERKSEAYLAS DLAERELTRD THKAIVFVK SVLLGGKSLK DLLPYGVIVC AFIFIPETAS VLDNVPVMIG NQKRPLTVAL IKYIAKSLNC DLVGDSYDTF YYCNSSAYGK NLISVSDNDF SNPQRALLSV GDLCYQAARS LHVAANYIR IFDRMPPGFQ PSKHLFRIIG VLDMETLTKM VTSNIAREPG MFCHDNVKDV LHRIGVYSPN HHFSAVILWK GWASTYAYMF NQEQLNMLSG TSGLAGDFGK YKLTYGSTFD EGVIVHVQYQF VTPEVVRKRN IYPDLSALKG GGS
<b>Source</b>	Yeast
<b>Protein Names</b>	Recommended name: Suppressor of RNA-mediated gene silencing Alternative name(s): Non-structural protein 10 Short name= Pns10
<b>Expression Region</b>	1-353
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.