



Recombinant Cyanidioschyzon merolae Pyruvate dehydrogenase E1 component subunit beta (pdhB)

Product Code	CSB-MP771221DZU
Abbreviation	pdhB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q85FX1
Product Type	Recombinant Protein
Immunogen Species	Cyanidioschyzon merolae (strain 10D) (Red alga)
Purity	≥85% (SDS-PAGE)
Sequence	MLHKLFMYEA LREAIDEEMA RDKRVFVLGE DVGHYGGSYK VTKQLHTKYG DLRVLDTPIA ENSFTGMAIG AAMTGLKPVV EGMNLSFLLL AFNQISNNAG MLHYTSGGNW SIPLVIRGPG GIGKQLSAEH SQRIEAYFQA VPGLKIVACS TPYNAKGLLK AAIRDNNPVL FLEHVLLYNL KQEIPKQEYV LPLDKAQVVR EGSDVTIITY SRMLHHVMQA VKQLVAQGMN PEVIDLISLK PIDLETLVTS VSKTHKAIIV EECMQTGGIA AEVMAQIYSH AFDELDAPIR RLSSKDVPTP YNGYLEQACL VQPTQIVEAV KTL MST
Source	Mammalian cell
Target Names	pdhB
Protein Names	Recommended name: Pyruvate dehydrogenase E1 component subunit beta EC= 1.2.4.1
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.