



# Recombinant Bovine Basigin (BSG)

<b>Product Code</b>	CSB-YP771275BO
<b>Abbreviation</b>	BSG
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q865R3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AAGTVFTTV EDLGSKILLT CSLNDSATEV TGHRLKGGV VLKEDALPGQ KTEFKVDSDD QWGEYSCVFL PEPMTANIQ LHGPPRVKAV KSSEHINEGE TAMLVCKSES VPPVTDWAWY KITDSEDKAL MNGSESRRFFV SSSQGRSELH IENLNMEADP GQYRCNGTSS KGSDQAIITL RVRSH
<b>Source</b>	Yeast
<b>Target Names</b>	BSG
<b>Protein Names</b>	Recommended name: Basigin Alternative name(s): EMMPRIN CD_antigen=CD147
<b>Expression Region</b>	22-205
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Extracellular domain
<b>Target Details</b>	This protein is a plasma membrane protein that is important in spermatogenesis, embryo implantation, neural network formation, and tumor progression. The encoded protein is also a member of the immunoglobulin superfamily. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>