



# Recombinant Lipoyl synthase, apicoplast (lipA)

<b>Product Code</b>	CSB-EP771386TOV-B
<b>Abbreviation</b>	lipA
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q86G50
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Toxoplasma gondii
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>GALSVSQ CSRRVSLGPL LSRVSSVSCT PSAAASALAS SLYPTDSLSS VEGSVAPRPP PSSLAFVLR VPPAAYSSSL SPSVLRFKHS LPRPLQGSV CAPGILGGAA GSARFAGCCG SQGRSCGSGK NPELPLKGSK DEVIPRVGTS TAGPRPDWFH VPAPQAASRG AEESRYQQLQ KQIRGLDLHT VCEEAKCPNI GECWNGGTAT LILLGDTCTR GCRFCAIKTS SKPPPPDPLE PEKVADAVAK WDIDYVVM TS VDRDDMPDGG AGHFARTVQL VKKAKPSMLI ECLVSDFQGM EESVRTLAQS GLDVYAHNIE TVRRLTPYVR DKRAKYDQSL RVLHLAKQFN PSLFTKSSIM LGLGETSEEV VRTLRLDRDH DVDVVTLGQY LRPTKQQLGV VEYVTPETFK KYQDIAEEMG FKYVASGPLV RSSYKAGEYY MKHLIDDARK HGRRETVKQV KLEADVGTKL GTTTTFQVNE KEA</p>
<b>Source</b>	E.coli
<b>Target Names</b>	lipA
<b>Protein Names</b>	<p>Recommended name: Lipoyl synthase, apicoplast EC= 2.8.1.8 Alternative name(s): Lipoate synthase Short name= LS Short name= Lip-syn Lipoic acid synthase</p>
<b>Expression Region</b>	64-543
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p>
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>