



# Recombinant Human RNA-binding Raly-like protein (RALYL)

<b>Product Code</b>	CSB-EP771422HU-B
<b>Abbreviation</b>	RALYL
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q86SE5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTGKTQTSNV TNKNDPKSIN SRVFIGNLNT AIVKKVDIEA IFSKYGKIVG CSVHKGYAFV QYMSERHARA AVAGENARVI AGQPLDINMA GEPKPYRPKP GNKRPLSALY RLESKEPFLS VGGYVFDYDY YRDDFYNRLF DYHGRVPPPP RAVIPLKRPR VAVTTTTRRGK GVFSMKGGSR STASGSTGSK LKSDELQTIK KELTQIKTKI DSSLGRLEKI EKQQKAEAEA QKKQLEESLV LIQEECVSEI ADHSTEEP AE GGP DADGEEM TDGIEEDFDE DGGHELFLQI K
<b>Source</b>	E.coli
<b>Target Names</b>	RALYL
<b>Protein Names</b>	Recommended name: RNA-binding Raly-like protein Short name= hRALYL Alternative name(s): Heterogeneous nuclear ribonucleoprotein C-like 3 Short name= hnRNP core protein C-like 3
<b>Expression Region</b>	1-291
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.