



Recombinant Human RNA-binding Raly-like protein (RALYL)

Product Code	CSB-MP771422HU
Abbreviation	RALYL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86SE5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTGKTQTSNV TNKNDPKSIN SRVFIGNLNT AIVKKVDIEA IFSKYGKIVG CSVHKGYAFV QYMSERHARA AVAGENARVI AGQPLDINMA GEPKPYRPKP GNKRPLSALY RLESKEPFLS VGGYVFDYDY YRDDFYNRLF DYHGRVPPPP RAVIPLKRPR VAVTTTTRRGK GVFSMKGGSR STASGSTGSK LKSDDELQTIK KELTQIKTKI DSSLGRLEKI EKQQKAEAEA QKKQLEESLV LIQEECVSEI ADHSTEEP AE GGP DADGEEM TDGIEEDFDE DGGHELFLQI K
Source	Mammalian cell
Target Names	RALYL
Protein Names	Recommended name: RNA-binding Raly-like protein Short name= hRALYL Alternative name(s): Heterogeneous nuclear ribonucleoprotein C-like 3 Short name= hnRNP core protein C-like 3
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.