



Recombinant Human Coiled-coil domain-containing protein 84 (CCDC84)

Product Code	CSB-MP771438HU
Abbreviation	CCDC84
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UT8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAPAQRCP LC RQTFFCGRGH VYSRKHQRQL KEALERLLPQ VEAARKAIRA AQVERYVPEH ERCCWCLCCG CEVREHLSHG NLTVLYGGLL EHLASPEHKK ATNKFWWENK AEVQMKEKFL VTPQDYARFK KSMVKGLDSY EEKEDKVIKE MAAQIREVEQ SRQEVVRSVL EPQAVPDPEE GSSAPRSWKG MNSQVASSLQ QPSNLDLPPA PELDWMETGP SLTFIGHQDI PGVGNIHSGA TPPWMIQDEE YIAGNQEIGP SYEEFLKEKE KQKLKLPDP RVGANFDHSS RTSAGWLPSF GRVWNNGRRW QSRHQFKTEA AAMKKQSHT E KS
Source	Mammalian cell
Target Names	CCDC84
Protein Names	Recommended name: Coiled-coil domain-containing protein 84
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.