



Recombinant Human Fermitin family homolog 3 (FERMT3)

Product Code	CSB-EP771441HU-B
Abbreviation	FERMT3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UX7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAGMKTASGD YIDSSWELRV FVGEEDPEAE SVTLRVGTGES HIGGVLLKIV EQINRKQDWS DHAIWWEQKR QWLLQTHWTL DKYGILADAR LFFGPQHRPV ILRLPNRRAL RLRASFSQPL FQAVAAICRL LSIRHPEELS LLRAPEKKEK KKKEKEPEEEE LYDLSKVVLA GGVAPALFRG MPAHFSDSAQ TEACYHMLSR PQPPDPPLL QRLPRPSSLS DKTQLHSRWL DSSRCLMQQG IKAGDALWLR FKYYSFFDL PKTDPVRLTQ LYEQARWDL LEEIDCTEEE MMVFAALQYH INKLSQSGEV GEPAGTDPGL DDL DVALSNL EVKLEGSAPT DVLDLSTTIP ELKDHLRIFR IPRRPRKLT L KGYRQHWVVF KETTL SYYKS QDEAPGDPIQ QLNLKGCEVV PDVNVSGQKF CIKLLVPSPE GMSEIYLRCQ DEQQYARWMA GCRLASKGRT MADSSYTSEV QAILAFLSLQ RTGSGGPGNH PHGPDASAEG LNPYGLVAPR FQRKFKAKQL TPRILEAHQN VAQLSLAEAEQ LRFIQAWQSL PDFGISYVMV RFKGSRKDEI LGIANNRLIR IDLAVGDVVK TWRFSNMRQW NVNWDIRQVA IEFDEHINVA FSCVSASCRI VHEYIGGYIF LSTRERARGE ELDEDLFLQL TGGHEAF
Source	E.coli
Target Names	FERMT3
Protein Names	Recommended name: Fermitin family homolog 3 Alternative name(s): Kindlin-3 MIG2-like protein Unc-112-related protein 2
Expression Region	1-667
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.