



Recombinant Human Serine/threonine-protein kinase pim-3 (PIM3)

Product Code	CSB-BP771444HU
Abbreviation	PIM3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86V86
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLLSKFGSLA HLCGPGGVDH LPVKILQPAK ADKESFEKAY QVGAVLGSSG FGTVYAGSRI ADGLPVAVKH VVKERVTEWG SLGGATVPLE VVLLRKVGAA GGARGVIRLL DWFERP DGFL LVLERPEPAQ DLFDFITERG ALDEPLARRF FAQVLA AVRH CHSCGVVHRD IKDENLLVDL RSGELKLIDF GSGALLKDTV YTFDGDTRVY SPPEWIRYHR YHGRSATVWS LGVLLYDMVC GDIPFEQDEE ILRGRLLFRR RVSPECQQLI RWCLSLRPSE RPSLDQIAAH PWMLGADGGV PESCDLRLCT LDPDDVASTT SSESSESL
Source	Baculovirus
Target Names	PIM3
Protein Names	Recommended name: Serine/threonine-protein kinase pim-3 EC= 2.7.11.1
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.