



Recombinant Human Uncharacterized protein C11orf74 (C11orf74)

Product Code	CSB-EP771446HU
Abbreviation	C11orf74
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86VG3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSAHMSGLEI MDEDQLIKDV LDKFLNCHEQ TYDEEFLNTF THLSQEDHVS KRGVFGTDSS ENIFTSAKVT HKNEADDYHL RNKTIFLRTS SQCLEEQVDN FLDLEDLDMD EEIKPQMSSED LLLLPGVEVEQ DVSTSIPSCI PFVAQPPTCE VKPKPSVKRM DKQTEEILGD EVQLFSLDEE FDYDNVMLTS KFSPAIEINI KELCKQQKRK DTSPDLEKSC D
Source	E.coli
Target Names	C11orf74
Protein Names	Recommended name: Uncharacterized protein C11orf74 Alternative name(s): Protein HEPIS
Expression Region	1-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.