



Recombinant Human Tax1-binding protein 1 (TAX1BP1)

Product Code	CSB-BP771448HU
Abbreviation	TAX1BP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86VP1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MTSFQEVLPLQ TSNFAHVIFQ NVAKSYLPNA HLECHYTLTP YIHHPKDWV GIFKVGWSTA RDYYTFLWSP MPEHYVEGST VNCVLAFAQGY YLPNDDGEFY QFCYVTHKGE IRGASTPFQF RASSPVEELL TMEDEGNSDM LVVTTKAGLL ELKIEKTMKE KEELLKLIIV LEKETAQLRE QVGRMERELN HEKERCDQLQ AEQKGLTEVT QSLKMENEEF KKRFSDATSK AHQLEEDIVS VTHKAIEKET ELDSLKDKLK KAQHEREQLE CQLKTEKDEK ELYKVHLKNT EIENTKLMSE VQTLKNLDGN KESVITHFKE EIGRLQLCLA EKENLQRTFL LTTSSKEDTC FLKEQLRKA E QVQATRQEV VFLAKELSDA VNVRDRTMAD LHTARLENEK VKKQLADAVA ELKLNAMKKD QDKTDTLEHE LRREVEDLKL RLQMAADHYK EKFKECQRLQ KQINKLSDQS ANNNNVFTKK TGNQQKV NDA SVNTDPATSA STVDVKPSPS AAEADFDIVT KGQVCEMTKE IADKTEKYNK CKQLLQDEKA KCNKYADELA KMELKWKEQV KIAENVKLEL AEVQDNYKEL KRSLENPAER KMEGQNSQSP QCFKTCSEQN GYVLTLSNAQ PVLQYGNPYA SQETRDGADG AFYPDEIQR PVRVPSWGLE DNVVCSQPAR NFSRPDGLD SEDSKEDENV PTAPDPPSQH LRGHGTGFCF DSSFDVHKKC PLCELMFPPN YDQSKFEEHV ESHWKVCPMC SEQFPDPYDQ QVFERHVQTH FDQNVLNFD</p>
Source	Baculovirus
Target Names	TAX1BP1
Protein Names	Recommended name: Tax1-binding protein 1 Alternative name(s): TRAF6-binding protein
Expression Region	1-789
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.