



# Recombinant Human Tax1-binding protein 1 (TAX1BP1)

<b>Product Code</b>	CSB-EP771448HU-B
<b>Abbreviation</b>	TAX1BP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q86VP1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MTSFQEVLPLQ TSNFAHVIFQ NVAKSYLPNA HLECHYTLTP YIHHPKDWV  GIFKVGWSTA RDYYTFLWSP MPEHYVEGST VNCVLAFAQGY YLPNDDGEFY  QFCYVTHKGE IRGASTPFQF RASSPVEELL TMEDEGNSDM LVVTTKAGLL  ELKIEKTMKE KEELLKLIIV LEKETAQLRE QVGRMERELN HEKERCDQLQ  AEQKGLTEVT QSLKMENEEF KKRFSDATSK AHQLEEDIVS VTHKAIEKET  ELDSLKDKLK KAQHEREQLE CQLKTEKDEK ELYKVHLKNT EIENTKLMSE  VQTLKNLDGN KESVITHFKE EIGRLQLCLA EKENLQRTFL LTTSSKEDTC  FLKEQLRKA E QVQATRQEV VFLAKELSDA VNVRDRTMAD LHTARLENEK  VKKQLADAVA ELKLNAMKKD QDKTDTLEHE LRREVEDLKL RLQMAADHYK  EKFKECQRLQ KQINKLSDQS ANNNNVFTKK TGNQQKV NDA SVNTDPATSA  STVDVKPSPS AA EADFDIVT KGQVCEMTKE IADKTEKYNK CKQLLQDEKA  KCNKYADELA KMELKWKEQV KIAENVKLEL AEVQDNYKEL KRSLENPAER  KMEGQNSQSP QCFKTCSEQN GYVLTLSNAQ PVLQYGNPYA  SQETRDGADG AFYPDEIQR PVRVPSWGLE DNVVCSQPAR NFSRPDGLD  SEDSKEDENV PTAPDPPSQH LRGHGTGFCF DSSFDVHKKC PLCELMFPPN  YDQSKFEEHV ESHWKVCPMC SEQFPPDYDQ QVFERHVQTH FDQNVLNFD</p>
<b>Source</b>	E.coli
<b>Target Names</b>	TAX1BP1
<b>Protein Names</b>	Recommended name: Tax1-binding protein 1 Alternative name(s): TRAF6-binding protein
<b>Expression Region</b>	1-789
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.