



Recombinant Human Caytaxin (ATCAY)

Product Code	CSB-EP771456HU-B
Abbreviation	ATCAY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86WG3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGTTTEATLRM ENVDVKEEWQ DEDLPRPLPE ETGVELLGSP VEDTSSPNT LNFNGAHRKR KTLVAPEINI SLDQSEGSLL SDDFLDTPDD LDINVDDIET PDETDSLEFL GNGNELEWED DTPVATAKMN PGDSADLFGD GTTEDGSAAN GRLWRTVIIG EQEHRIDLHM IRPYMKVVTH GGYGEGEGLNA IIVFAACFLP DSSLPDYHYI MENLFLYVIS SLELLVAEDY MIVYLNATP RRRMPGIGWL KKCYQMIDRR LRKNLKSLLI VHPSWFIRTV LAISRPFISV KFINKIQYVH SLEDLEQLIP MEHVQIPDCV LQYEEERLKA RRESARPQPE FVLPRSEEKP EVAPVENRSA LVSEDQETSM S
Source	E.coli
Target Names	ATCAY
Protein Names	Recommended name: Caytaxin Alternative name(s): Ataxia cayman type protein BNIP-H
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a neuron-restricted protein that contains a CRAL-TRIO motif common to proteins that bind small lipophilic molecules. Mutations in this gene are associated with cerebellar ataxia, Cayman type.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.