



Recombinant Human Uncharacterized protein C19orf68 (C19orf68)

Product Code	CSB-MP771465HU
Abbreviation	C19orf68
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86XI8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MERPEPPPGT AAGQEEQELR ERAFFSWAEF SRFFDAWCQQ RLALFFVKSS MHLARCRWAS APPLYTLIDV LKYSYVRLVC KDVRAPSRPA VGPPQPGCPA FIIVKLSPLR DRLVVTECQL THSHPACPLE FAYYFRPGHL LANACLVRT TNKISKQFVA PADVRRLLSY CKGRDHGULD ALHVLEGLFR TDPEAKVKLV FVEDQAVVET VFFLTSRTRA LLRRFPRMLL VDRLPGLQGA LDLLAVLCVD GSGRARQAAC CVARPGTPSL LRFALASLLQ SAPDVKGRVR CLTAGPEVAA QLPAVRQLLP CARVQICRAQ GLETLFKAQ ELGGAGREDP GLWSRLCRLA GASSPAAYDE ALAELHAHGP AAFVDYFERN WEPRRDMWVR FRAFEAARDL DACALVRGHR RLLRRLSPS RGVAQCLRDV VAMQWADAAG EAVPEGPDGG GPWLEDEPGR GAQGENERVR GLETGDWGGA PKEGSIWRGA QMEKEWARAL ETRDWGGAQF EGEKGRALQI RDWRGGRLN QKPRGLEGGV LRGSKLEKLG LRGPEIRDWR GPQLEGEKDW GLEGYVWRGS QLEDQALRGL EGYTWRVAQL EDRRSTDLR GTQFDYER GK GESTEDR
Source	Mammalian cell
Target Names	ZSWIM9
Protein Names	Recommended name: Uncharacterized protein C19orf68
Expression Region	1-627
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.