



# Recombinant *Oryza sativa* subsp. japonica Protochlorophyllide reductase A, chloroplastic (PORA)

<b>Product Code</b>	CSB-EP771939OFG-B
<b>Abbreviation</b>	PORA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7XKF3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QKAKK PSLVVRAVAT PAAPVASPGA GTSKADGKKT LRQGVVITG ASSGLGLAAA KALAETGKWH VVMACRDFLK AATAAKAAGM AAGSYTMHL DLASLDSVRQ FVDNFRMSGM PLDALVCNAA IYRPTARQPT FTADGYEMSV GVNHLGHFLL ARLMLDDLKK SDYPSRRLII LGSITGNTNT LAGNVPPKAG LGDLRGLAGG LRGQNGSAMI DGAESFDGAK AYKDSKICNM LTMQEFHRRF HEETGITFAS LYPGCIATTG LFREHIPLFR LLFPPFQRFV TKGFVSEAES GKRLAQVVG D PSLTKSGVYW SWNKDSASFE NQLSQEASDP EKARKLWDLS EKLVLV
<b>Source</b>	E.coli
<b>Target Names</b>	PORA
<b>Protein Names</b>	Recommended name: Protochlorophyllide reductase A, chloroplastic Short name= PCR A EC= 1.3.1.33 Alternative name(s): NADPH-protochlorophyllide oxidoreductase A Short name= POR A
<b>Expression Region</b>	36-387
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.