



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 30 (Os04g0663200, LOC_Os04g56750)

Product Code	CSB-EP771941OFG-B
Abbreviation	Os04g0663200, LOC_Os04g56750
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XM16
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MMGSRRSRRV SWASGGNLCK VRLFLSEDSP SQAGLRPQDN LQAKGSWLLH AAGPSSDDSL PPGFESLPPS NDLKIDMSQI PLIRWKCP IVLEQDWHIV AGEESREIEI QNERINGALE AIYPRPSNIP PNPFLSLDVK DAHYYDDSKTL LVPLIPLEDD DASDQLEGPT LDLPSHYNIT GVSNTPVSAE QQPPCGGAIS SGFTIEPQAA VSATVTAIMQ TIQSNQNGSM ADQNGSTIDQ ELLFKILSDP SQLQRLMKEC GPVRHEQSAS SSVVAPLVSI PPPQITASSP APFSDHVGTF HGMNPTLPPP PPMNRPST IPSVAMNHPP SSSPAMNFGS ALPSSSPSVN FGSVPGRGVG YYKTLIHQHG GERLEQPF EQ HGMQFGMYRQ PGPPQNGGID AMNGAAAMVS RDGKVRPMKP CAYFNPKGC RENGASCTFLH DASAPTRKDH QKQKGSKRIK LDNTMGGRN
Source	E.coli
Target Names	Os04g0663200
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 30 Short name= OsC3H30
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.