



Recombinant *Oryza sativa* subsp. japonica Anamorsin homolog 1 (Os04g0674400, LOC_Os04g57810)

Product Code	CSB-MP771945OFG
Abbreviation	Os04g0674400, LOC_Os04g57810
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XQ97
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAATVAEALA VTDELALPLR AVGDIAAAAG VSREEVVIT QCASLGGKLP FDDASVGSVL AVIKKVENLG DLFITEISR V LKAGGMVLIQ SSPSDQDPNN SIQRKLLGG FVDVQASAAS SQDSEHSVTI KAKKVSWSLG SSFPLKKATK GLPKIQIDDD SELIDEDSLL TEDDLKKPEL PVVGDCEVGA TRKACKNCTC GRAEAEEKVE KLNLTSEQIN NPQSACGNCG LGDAFRCGTC PYRGLPAFKP GEKIALPGNF LAADM
Source	Mammalian cell
Target Names	Os04g0674400
Protein Names	Recommended name: Anamorsin homolog 1 Alternative name(s): Fe-S cluster assembly protein DRE2 homolog 1
Expression Region	1-265
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.