



Recombinant *Oryza sativa* subsp. japonica PHD finger protein ALFIN-LIKE 4 (Os04g0444900, LOC_Os04g36730)

Product Code	CSB-EP771955OFG-B
Abbreviation	Os04g0444900, LOC_Os04g36730
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XUW3
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDGGYGSVTI VHDARSPEDV FQDFCGRRSG IVKALTIEVE KFYKQCDPEK ENLCLYGLPN GTWAVTLPAD EVPPELPEPA LGINFARDGM QEKDWLSLIA VHSDSWLLSV AFYFGARFGF DKKARERLFM MTSSLPTVFE VVSGGVNTQS KTANGSSKNK SGSKPPKRPN SDSKPQKQVQ AKYEEENGGR GNGGDEDQAE TICGACGEAY ANGEFWICCD ICETWFHGKC VRITPAKAEH IKHYKCPGCS NKRTRE
Source	E.coli
Target Names	Os04g0444900
Protein Names	Recommended name: PHD finger protein ALFIN-LIKE 4
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.