



Recombinant Enterobacteria phage RB69 Deoxynucleotide monophosphate kinase (1)

Product Code	CSB-BP771963ELM
Abbreviation	1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Y4Y9
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage RB69 (Bacteriophage RB69)
Purity	>85% (SDS-PAGE)
Sequence	MELIFLSGIK RSGKDTTADY INSNFKSIKY QLAYPIKDAL AIAWERKHAE NPDVFTELKY EYFEGIGYDR ETPLNLNKLD VIELMEETLI YLQRQYLPIN GVNILSSLEG GYSYLDIKPY EALREAINNI NDTWSIRRLM QALGTDVVVN LFDRMYWVKL FALNYMDYIG SDFDYVVTD TRQVHEMETA RAMGATVIHV VRSGTESTDK HITEAGLPIE EGDLVITNDG SLEELYSKIE KILR
Source	Baculovirus
Target Names	1
Protein Names	Recommended name: Deoxynucleotide monophosphate kinase Short name= DNK Short name= dNMP kinase EC= 2.7.4.13 Alternative name(s): Gp1
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.