



Recombinant Human L-xylulose reductase (DCXR)

Product Code	CSB-MP772022HU
Abbreviation	DCXR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z4W1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MELFLAGRRV LVTGAGKGIG RGTVQALHAT GARVVAVSRT QADLDSLVR CPGIEPVCVD LGDWEATERA LGSVGPVDLL VNNAVALLQ PFLEVTKEAF DRSFEVNLRA VIQVSQIVAR GLIARGVPGA IVNVSSQCSQ RAVTNHSVYC STKGALDMLT KVMALGLPH KIRVNAVNP VVMTSMGQAT WSDPHKAKTM LNRIPLGKFA EVEHVVNAIL FLLSDRSGMT TGSTLPVEGG FWAC
Source	Mammalian cell
Target Names	DCXR
Protein Names	Recommended name: L-xylulose reductase Short name= XR EC= 1.1.1.10 Alternative name(s): Carbonyl reductase II Dicarbonyl/L-xylulose reductase Kidney dicarbonyl reductase Short name= kiDCR Sperm surface protein P34H
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.