



Recombinant Human Abhydrolase domain-containing protein 12B (ABHD12B)

Product Code	CSB-EP772026HU-B
Abbreviation	ABHD12B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z5M8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDAQDCQAAA SPEPPGPPAR SCVAAWDMV DRNLRYFPHS CSMLGRKIAA LYDSFTSKSL KEHVFLPLID MLIYFNFFKA PFLVDLKKPE LKIPHTVNFY LRVEPGVMLG IWHTVPSCRG EDAKGDCCW YEAALRDGNP IIVYLHGSAE HRAASHRLKL VKVLSGGGFH VLSVDYRGFG DSTGKPTEEG LTTDAICVYE WTKARSGITP VCLWGHSLGT GVATNAAKVL EEKGCVPDAI VLEAPFTNMW VASINYPLLK IYRNIPGFLR TLMDALRKDK IIFPNDEVK FLSSPLLILH GEDDRTVPLE YGKKLYEIAR NAYRNKERVK MVIFPPGFQH NLLCKSPTLL ITVRDFLSKQ WS
Source	E.coli
Target Names	ABHD12B
Protein Names	Recommended name: Abhydrolase domain-containing protein 12B EC= 3.-.-.
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.