



# Recombinant Danio rerio Keratin, type I cytoskeletal 18 (krt18)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP772048DIL-B   |
| <b>Abbreviation</b>      | krt18   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q7ZTS4  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Danio rerio (Zebrafish) (Brachydanio rerio)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | SLRTSYSVR SSTSQVPVSQ VSIKRTTNVP TYRAASIYGG AGGQGTRISS<br>ASYSGVRSGL GVPSMSSSIQ VSASGSTGEI MGNEKMAMQN LNDRLASYLE<br>KVRILEQANS KLELKIREAL EKRGPDVHDY SRFQPIVDEL RKKIFDATTN<br>NARLVLQIDN ARLAADDFRV KYEELSIRQ GVEADITGLR KVIDDTNLNR<br>MNLESEIEAL KEELIFLKKH HDNEVMELRN QISQSGVQVD VDAPKGQDLS<br>QIMEEIRAKY EKMALKNQEE LKAWHESQIT EVQVQVTQNT EALQGARSEV<br>NELRRQIQL EIELESQKNL KGSLEGTLRD TEMRYNMEIE NLNTIILQLE<br>AELTQLRGNI QHQTQEYEAL LNIKMKLEAE IATYRRLLDG GDFKLQDALE<br>EQKKVKVMTV TQTLVDGKVV SSSTETKERK L |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | krt18   |
| <b>Protein Names</b>     | Recommended name: Keratin, type I cytoskeletal 18 Alternative name(s):<br>Cytokeratin-18 Short name= CK-18 Dorsal aorta proneprin kinesin-1 DreK18<br>Keratin-18 Short name= K18  |
| <b>Expression Region</b> | 2-431   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.