



Recombinant Danio rerio Protein rogdi homolog (rogdi)

Product Code	CSB-BP772061DIL
Abbreviation	rogdi
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZVT5
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MTAASQAERT VLEEEFNWLL KEEVHAVLKQ LQDILKEASR RLSMPSPGLE GQLKQENFIL GSSTMDQVKG VLTLQGEALT QADINIKVAK SSQVMHFAFR DDKQWKLQIQI QDARNHVNQA LQLLSSRDDS YHFKTGAEVN KLMDAIMLQL TRARNRLTTP ASMTLPELAA SGLMKMFTPP MPGDVMVNFY INLSKLCLTV YQLHVLQPNT TKNFKPAGSS VLHNPGAMFE HNNTKFEVSH VHKVECVVPW LNDTLVFFTI SLQLCQQLKD KISVFSSFWN YRPF
Source	Baculovirus
Target Names	rogdi
Protein Names	Recommended name: Protein rogdi homolog
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.