



Recombinant Danio rerio Protein FAM3C (fam3c)

Product Code	CSB-BP772085DIL
Abbreviation	fam3c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZYY4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	NMNFNLGKVF TRYAPDAVP TRPLRHKCGL SKSCPEDHFA FKITSGAASV VGPKMCFQDN VIMSGVKNNI GRGINIALLN GKTGELTKTD SFDMWTGDVN LLIKFLKDIE DGSIVMMATF DDPATKMNDE ARNLIADLGS SSIILGFRD NWVFGGKGI KTKSPFEQHI KNNAAETNKYE GWPEVLEMEG CIPIKHD
Source	Baculovirus
Target Names	fam3c
Protein Names	Recommended name: Protein FAM3C
Expression Region	31-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.