



Recombinant Danio rerio O-phosphoseryl-tRNA (Sec) selenium transferase (sepsecs)

Product Code	CSB-MP772095DIL
Abbreviation	sepsecs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q803A7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MNSENFRLLSE KLVSASYVRQ GLDARRGHEQ LIRRLLDQGK CPEEGWSEST IELFLNELAV MDSNNFLGNC GVGEREGRVA SSLVARRHYR LIHGIGRSGD IAAVQPKAAG SLLNKITNS VVLDVLKLTG VRSVSSCFVV PMATGMSLTL CFLTLRHRRP KARYILWPRI DQKSCFKSMV TAGFEPVVIE NVLEGDELRT NLEEVERKIE EFGAENTLCV HSTTSCFAPR VPDRLEELSV LCAKHDIPHI VNNAYGVQPS KCMHLIQQGA RVGRIDAFVQ SLDKNFVVPV GGAIAGFDE NFIQEISKIY PGRASASPSL DVLITLLTLG ANGYKLLAD RKELYGHFAQ ELSALAARHG ERLKTPHNP ISLAMSLNHL EAHSSSAVTQ LGSMLFTRQV SGARVVPLGV QQTVSGHTFS GFMSHSEAYP CPYLNAASAI GITKGDVTVK MKRLGKCLKI LKKEKSDPAD LEAEDGELEE SPQRSTELRV
Source	Mammalian cell
Target Names	sepsecs
Protein Names	Recommended name: O-phosphoseryl-tRNA(Sec) selenium transferase EC=2.9.1.2 Alternative name(s): Selenocysteine synthase Short name= Sec synthase Selenocysteinyl-tRNA(Sec) synthase Sep-tRNA:Sec-tRNA synthase Short name=
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.