



# Recombinant Xaa-Pro dipeptidase (pepQ)

<b>Product Code</b>	CSB-BP772757SZB
<b>Abbreviation</b>	pepQ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q83PG0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Shigella flexneri
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESLASLYKN HIATLQERTR DALARFKLDA LLIHSGELFN VFLDDHPYPF KVNPQFKAWV PVTQVPNCWL LVDGVNPKPL WFYLPVDYWH NVEPLPTSFV TEDVEVIALP KADGIGSLLP AARGNIGYIG PVPERALQLG IEASNINPKG VIDYLHYYS FKTEYELACM REAQKMAVNG HRAAEEAFRS GMSEFDINIA YLTATGHRDT DVPYSNIVAL NEHA AVLHYT KLDHQAPEEM RSFLLDAGAE YNGYAADLTR TWSAKSDNDY AQLVKDVNDE QLALIATMKA GVS YVDYHIQ FHQRIAKLLR KHQITDMSE EAMVENDLTG PFMPHGIGHP LGLQVHDVAG FMQDDSGTHL AAPAKYPYLR YTRILQPGMV LTIEPGIYFI ESLLAPWREG QFSKHFNWQK IEALKPFGGI RIEDNVVIHE NNVENMTRDL KLA
<b>Source</b>	Baculovirus
<b>Target Names</b>	pepQ
<b>Protein Names</b>	Recommended name: Xaa-Pro dipeptidase Short name= X-Pro dipeptidase EC= 3.4.13.9 Alternative name(s): Imidodipeptidase Proline dipeptidase Short name= Prolidase
<b>Expression Region</b>	1-443
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.