



Recombinant *Oryza sativa* subsp. *japonica* 12-oxophytodienoate reductase 1 (OPR1)

Product Code	CSB-MP772830OFG
Abbreviation	OPR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84QK0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MVHAPAKVAA AAAILLTPY KMGQLELSHR VVLAPLTRCR SYGNVPQPHA AVVYSQRATR GLLIAEATD ISPTAQGYPE TPGIYTQQQI EAWKPIVDAV HRKGALFFLQ IWHVGRVSTT DFQPNGQAPI SSTDKQITPD DSGMVYSKPR RLRTDEIPQI IDDFRRAARN AIEAGFDGVE IHGAHGYLE QFMKDSANDR TDEYGGLEN RCRFAVEVID AVVAEVGAHR VGIRLSPFVD FMDCFDSDPV ALGSYMQQL NKHPGFLYCH MVEPRMAIIE GRRKIAHGILL PFRKQFNGTF IAAGGYDREE GNKVVADGYA DLVAYGRLFL ANPDLPRRFE LDAPLNRYDR STFYTQDPVV GYTDYPFLEE IDEESRTTYA
Source	Mammalian cell
Target Names	OPR1
Protein Names	Recommended name: 12-oxophytodienoate reductase 1 EC= 1.3.1.42 Alternative name(s): 12-oxophytodienoate-10,11-reductase 1 Short name= OPDA-reductase 1 Short name= OsOPR1
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.