



Recombinant *Oryza sativa* subsp. japonica Zinc finger A20 and AN1 domain-containing stress-associated protein 7 (SAP7)

Product Code	CSB-MP772856OFG
Abbreviation	SAP7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q852K6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MASMKRKCPD DETACGSGAG AAMCVTGCGF FGSEATNNMC SRCYREHSAD NDAVEEAAAA NSDLELVGVA ETTTKKARMS AVVPVAVASS SSAAAEQPAA KAATAPNRCA ACRKKVGLTG FKCRCGGNFC GGHRHADAHG CGFDYKSAGK EQIAKQNPLV VADKLATRI
Source	Mammalian cell
Target Names	SAP7
Protein Names	Recommended name: Zinc finger A20 and AN1 domain-containing stress-associated protein 7 Short name= OsSAP7
Expression Region	1-169
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.