



# Recombinant Adiantum capillus-veneris Maturase K (matK)

<b>Product Code</b>	CSB-MP772884AXO
<b>Abbreviation</b>	matK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q85FP0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Adiantum capillus-veneris (Maidenhair fern)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKVNYGSPAK SGPLHKNGEF YINKNCFLHP LLFLSEENFY LTDGKRRSHG ADANLVFGAW STVAVKRLIG SVRDPNLNFL KIYDSVSVRN LTDGLDVDLY LHPLLKMTYL ILGIALFPKI RAETSSKSKM LQSIHSMFLF LEDRFSKSNH ILEADLPHNL HLETLIRLFR RQIKDVSFLH LLRIVFRKRK IFCGKTFYSP GGGQDGSVDI PVRNFYIFEI DSLLLIPWKQ VYKFRVNYLS PIDSCNIIRK EIYASAYKFK WNKASIDYSF SRSLWIHYGR WRNKFLIASE GTHYFVKKML YYLWILLKYH FHYRIKSNEP WIRKLLPTSC VSFLGYTLA QLVSKNVRIE TVTDLYISIL GGKGFYPKIP NSIIITTLAK QRFCDFTRP IGKSAWVTST DDKIIDGYVQ LWQVFSLYYG ASMNQYRLRR LIFLLQMSCD STLAGKHRST IRLLRCKSNV EALNQILASR KFELSSSRRV WRSSSIRSVL VQFTVLDIGL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	matK
<b>Protein Names</b>	Recommended name: Maturase K Alternative name(s): Intron maturase
<b>Expression Region</b>	1-500
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.