



# Recombinant Human Serine/threonine-protein kinase 32C (STK32C)

<b>Product Code</b>	CSB-EP773035HU
<b>Abbreviation</b>	STK32C
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q86UX6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MRSGAERRGS SAAASPGSPP PGRARPAGSD APSALPPPAA GQPRARDSGD VRSQPRPLFQ WSKWKKRMGS SMSAATARRP VFDDKEDVNF DHFQILRAIG KGSFGKVCIV QKRDTEKMYA MKYMNKQQCI ERDEVNRNFR ELEILQEIEH VFLVNLWYSF QDEEDMFMVV DLLLGGDLRY HLQQNVQFSE DTVRLYICEM ALALDYLRGQ HIIHRDVKPD NILLDERGHA HLTDFNIATI IKDGERATAL AGTKPYMAPE IFHSFVNGGT GYSFEVDWWS VGVMAYELLR GWRPYDIHSS NAVESLVQLF STVSVQYVPT WSKEMVALLR KLLTVNPEHR LSSLQDVQAA PALAGVLWDH LSEKRVEPGF VPKNGRLHCD PTFELEEMIL ESRPLHKKKK RLAKNKS RDN SRDSSQSEND YLQDCLDAIQ QDFVIFNREK LKRSQDLPRE PLPAPESRDA AEPVEDEAER SALPMC GPIC PSAGSG
<b>Source</b>	E.coli
<b>Target Names</b>	STK32C
<b>Protein Names</b>	Recommended name: Serine/threonine-protein kinase 32C Short name= PKE EC= 2.7.11.1 Alternative name(s): YANK3
<b>Expression Region</b>	1-486
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is a member of the serine/threonine protein kinase family. The specific function of this kinase is not known.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.