



Recombinant Human Serine/threonine-protein kinase 32C (STK32C)

Product Code	CSB-MP773035HU
Abbreviation	STK32C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UX6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRSGAERRGS SAAASPGSPP PGRARPAGSD APSALPPPAA GQPRARDSGD VRSQPRPLFQ WSKWKKRMGS SMSAATARRP VFDDKEDVNF DHFQILRAIG KGSFGKVCIV QKRDTEKMYA MKYMNKQQCI ERDEVRNVFR ELEILQEIEH VFLVNLWYSF QDEEDMFMVV DLLLGGDLRY HLQQNVQFSE DTVRLYICEM ALALDYLRGQ HIIHRDVKPD NILLDERGHA HLTDFNIATI IKDGERATAL AGTKPYMAPE IFHSFVNGGT GYSFEVDWWS VGV MAYELLR GWRPYDIHSS NAVESLVQLF STVSVQYVPT WSKEMVALLR KLLTVNPEHR LSSLQDVQAA PALAGVLWDH LSEKRVEPGF VPNKGR LHCD PTFELEEMIL ESRPLHKKKK RLAKNKS RDN SRDSSQSEND YLQDCLDAIQ QDFVIFNREK LKRSQDLPRE PLPAPESRDA AEPVEDEAER SALPMCGPIC PSAGSG
Source	Mammalian cell
Target Names	STK32C
Protein Names	Recommended name: Serine/threonine-protein kinase 32C Short name= PKE EC= 2.7.11.1 Alternative name(s): YANK3
Expression Region	1-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the serine/threonine protein kinase family. The specific function of this kinase is not known.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.