



# Recombinant Human Protein FAM160B2 (FAM160B2)

<b>Product Code</b>	CSB-BP773037HU
<b>Abbreviation</b>	FAM160B2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q86V87
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>           MLSRLGALLQ EAVGAREPSI DLLQAFVEHW KGITHYYIES TDESTPAKKT            DIPWRLKQML DILVYEEQQQ AAAGEAGPCL EYLLQHKILE TLCTLGKAEY            PPGMRQQVFQ FFSKVLAQVQ HPLLHYLSVH RPVQKLLRLG GTASGSVTEK            EEVQFTTVLC SKIQQDPELL AYILEGKKIV GRKKACGEPT ALPKDTTSHG            DKDCSHDGAP ARPQLDGESC GAQALNSHMP AETEELDGGT TESNLITSL            GLCQSKKSRV ALKAQENLLL LVSMASPAAL TYLVQSSACC PAIVRHLCQL            YRSMPVFLDP ADIATLEGIS WRLPSAPSDE ASFPGKEALA AFLGWFDYCD            HLITEAHTVV ADALAKAVAE NFFVETLQPQ LLHVSEQSIL TSTALLTAML            RQLRSPALLR EAVAFLLGTD RQPEAPGDNP HTLYAHLIGH CDHLSDEISI            TTLRLFEEEL QKPHEGIIHS LVLRNLEGRP YVAWGSPEPE SYEDTLDLEE            DPYFTDSFLD SGFQTPAKPR LAPATSYDGK TAVTEIVNSF LCLVPEEAKT            SAFLEETGYD TYVHDAYGLF QECSSRVASW GWPLTPTPLD PHEPERPFFE            GHFLRVLFDR MSRILDQPYS LNLQVTSVLS RLALFPHPHI HEYLLDPYIS            LAPGCRSLFS VLVRVIGDLM QRIQRVPQFP GKLLLVKRLQ TGQAPGEQLD            HQTLLQGVVV LEEFCKELAA IAFVKFPPHD PRQNVSPAPE GQV         </p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	FAM160B2
<b>Protein Names</b>	Recommended name: Protein FAM160B2 Alternative name(s): Retinoic acid-induced protein 16
<b>Expression Region</b>	1-743
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.