



Recombinant Human Protein FAM160B2 (FAM160B2)

Product Code	CSB-MP773037HU
Abbreviation	FAM160B2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86V87
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MLSRLGALLQ EAVGAREPSI DLLQAFVEHW KGITHYYIES TDESTPAKKT DIPWRLKQML DILVYEEQQQ AAAGEAGPCL EYLLQHKILE TLCTLGKAEY PPGMRQQVFQ FFSKVLAQVQ HPLLHYSVH RPVQKLLRLG GTASGSVTEK EEVQFTTVLC SKIQQDPELL AYILEGKKIV GRKKACGEPT ALPKDTSHG DKDCSHDGAP ARPQLDGESC GAQALNSHMP AETEELDGGT TESNLITSL GLCQSKKSRV ALKAQENLLL LVSMASPAAL TYLVQSSACC PAIVRHLCQL YRSMPVFLDP ADIATLEGIS WRLPSAPSDE ASFPGKEALA AFLGWFDYCD HLITEAHTVV ADALAKAVAE NFFVETLQPQ LLHVSEQSIL TSTALLTAML RQLRSPALLR EAVAFLLGTD RQPEAPGDNP HTLYAHLIGH CDHLSDEISI TTLRLFEEEL QKPHEGIIHS LVLRNLEGRP YVAWGSPEPE SYEDTLDLEE DPYFTDSFLD SGFQTPAKPR LAPATSYDGK TAVTEIVNSF LCLVPEEAKT SAFLEETGYD TYVHDAYGLF QECSSRVASW GWPLTPTPLD PHEPERPFFE GHFLRVLFDR MSRILDQPYS LNLQVTSVLS RLALFPHPHI HEYLLDPYIS LAPGCRSLFS VLVRVIGDLM QRIQRVPQFP GKLLLVKRLQ TGQAPGEQLD HQTLLQGVVV LEEFCKELAA IAFVKFPPHD PRQNVSPAPE GQV </p>
Source	Mammalian cell
Target Names	FAM160B2
Protein Names	Recommended name: Protein FAM160B2 Alternative name(s): Retinoic acid-induced protein 16
Expression Region	1-743
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.