



Recombinant Human Protein FAM90A1

(FAM90A1)

Product Code	CSB-MP773059HU
Abbreviation	FAM90A1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86YD7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMARRDPKPG AKRLVRAQTL QKQRRAPVGP RAPPPDEEDP RLKCKNCEAF GHTARSTRCP MKCWKAALVP PNFGEKEGKE NLKPWKQVE ANPGPLNKDK GEKEERPRPQ DPQRKALLHI FSRKPPEKPL PNQKGSTESS DYLRVASGPM PVHTTSKRPR VDPVLSDRSA TEMSDRGSVL ASLSPLRKAS LSSSSSLGPK ERQTGAAADI PQTAVRHQGP EPLLVKPTH SSPAGGCREV PQAASKTHGL LQAVSPQAQD KRPAVTSQPC PPAATHSLGL GSNLSFGPGA KRSAPAPIQA CLNFPKKPRL GPFQIPESAI QGGELGAPEN LQPPPAATEL GPRTSPQTGT RTPAQVLSGD RQPPHSRPCL PTAQACTMSH HPAASHDGAQ PLRVLFRRLE NGRWSSSLT APSFHSPEKP GAFLAQSPHV SEKSEGPCVR VPPSVLYEDL QVPSSSEDSD SDLE
Source	Mammalian cell
Target Names	FAM90A1
Protein Names	Recommended name: Protein FAM90A1
Expression Region	1-464
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.