



Recombinant Human Decaprenyl-diphosphate synthase subunit 2 (PDSS2)

Product Code	CSB-BP773060HU
Abbreviation	PDSS2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86YH6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNFRQLLLHL PRYLGASGSP RRLWWSPSLD TISSVGSWRG RSSKSPAHWN QVVSEAEKIV GYPTSFMSLR CLLSDELSNI AMQVRKLVGT QHPLTTARG LVHDSWNSLQ LRGLVLLIS KAAGPSSVNT SCQNYDMVSG IYSCQRSLAE ITELIHIAL VHRGIVNLNE LQSSDGPLKD MQFGNKIAIL SGDFLLANAC NGLALLQNTK VVELLASALM DLVQGVYHEN STSKESYITD DIGISTWKEQ TFLSHGALLA KSCQAAMELA KHDAEVQNMA FQYGKHMAMS HKINSVDVQPF IKEKTSDSMT FNLNSAPVVL HQEFLGRDLW IKQIGEAQEK GRLDYAKLRE RIKAGKGVTS AIDLCRYHGN KALEALESFP PSEARSALEN IVFAVTRFS
Source	Baculovirus
Target Names	PDSS2
Protein Names	Recommended name: Decaprenyl-diphosphate synthase subunit 2 EC=2.5.1.91 Alternative name(s): All-trans-decaprenyl-diphosphate synthase subunit 2 Candidate tumor suppressor protein Decaprenyl pyrophosphate synthase subunit 2
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.