



# Recombinant Human Decaprenyl-diphosphate synthase subunit 2 (PDSS2)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP773060HU                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Abbreviation</b>      | PDSS2                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                          |
| <b>Uniprot No.</b>       | Q86YH6                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MNFRQLLLHL PRYLGASGSP RRLWWSPSLD TISSVGSWRG<br>RSSKSPAOWN QVVSEAEKIV GYPTSFMSLR CLLSDELSNI AMQVRKLVGT<br>QHPLTTARG LVHDSWNSLQ LRGLVLLIS KAAGPSSVNT SCQNYDMVSG<br>IYSCQRLAE ITELIHIAL VHRGIVNLNE LQSSDGPLKD MQFGNKIAIL<br>SGDFLLANAC NGLALLQNTK VVELLASALM DLVQGVYHEN STSKESYITD<br>DIGISTWKEQ TFLSHGALLA KSCQAAMELA KHDAEVQNMA FQYGKHMAMS<br>HKINSVDVQPF IKEKTSDSMT FNLNSAPVVL HQEFLGRDLW IKQIGEAQEK<br>GRLDYAKLRE RIKAGKGVTS AIDLCRYHGN KALEALESFP PSEARSALEN<br>IVFAVTRFS |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Target Names</b>      | PDSS2                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: Decaprenyl-diphosphate synthase subunit 2 EC=2.5.1.91 Alternative name(s): All-trans-decaprenyl-diphosphate synthase subunit 2 Candidate tumor suppressor protein Decaprenyl pyrophosphate synthase subunit 2                                                                                                                                                                                                                                             |
| <b>Expression Region</b> | 1-399                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                 |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.                                                                                                                                                                                                                                                                                                                                  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.