



# Recombinant Arabidopsis thaliana Probable protein phosphatase 2C 35 (At3g06270)

<b>Product Code</b>	CSB-MP773514DOA
<b>Abbreviation</b>	At3g06270
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7XJ53
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGCVQCKCCS RYPSSSSDGD SRGPLEANGV LKGKDQKPLG SIHVPSPNFD MVYSVLSQRG YYPDSPDKEN QDTYCIKTEL QGNPNVHFFG VFDGHGVLGT QCSNFVKERV VEMLSEDPTL LEDPEKAYKS AFLRVNEELH DSEIDDSMSG TTAITVLVVG DKIVVANVGD SRAVLAVKDR NRILAEDLSY DQTPFRKDEC ERVKACGARV LSV DQVEGLK DPNIQTWANE ESEGGDPPRL WVQNGMYPGT AFTRSVGDFT AESIGVIAEP EVSMVHLSPN HFFFVVASDG IFEFLPSQAV VDMVGRYADP RDGCAAAA AE SYKLWLEHEN RTDDITIIIV QIKKLSNE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At3g06270
<b>Protein Names</b>	Recommended name: Probable protein phosphatase 2C 35 Short name= AtPP2C35 EC= 3.1.3.16
<b>Expression Region</b>	1-348
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.