



Recombinant *Oryza sativa* subsp. *japonica* DEAD-box ATP-dependent RNA helicase 6 (Os04g0533000, LOC_Os04g45040)

Product Code	CSB-EP773520OFG-B
Abbreviation	Os04g0533000, LOC_Os04g45040
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XMK8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDPRARYPPG IGNGRGGNPN YYGRGPPPSQ HQQHQHQHQQQ PPHPHHHQYV QRQPQPQQTP HNSQHQQWLR RNQIAAEAAG ASEQKAPPVA DGIDSSSQDW KAQLKLPPQD TRYRTEDVTA TKGNEFEDYF LKRELLMGIY EKGFERPSPI QEESIPIALT GSDILARAKN GTGKTAAFCI PALEKIDQDK NAIQVVILVP TRELALQTSQ VCKELGKHLK IQVMVTTGGT SLKDDIVRLY QPVHLLVGTP GRILDLTKKG VCVLKNCSML VMDEADKLLS PEFQPSIQEL IRYLPSNRQI LMFSATFPVT VKEFKDKYLP KPYVINLMDE LTLKGITQFY AFVEERQKVH CLNTLFSKLQ INQSIIFCNS VNRVELLAKK ITELGYSCFY IHAKMLQDHR NRVFHDFRNG ACRNLVCTDL FTRGIDIQAV NVVINDFDPK SAETYLHRVG RSGRFGHLGL AVNLITYEDR FNLYRIEQEL GTEIKPIPPQ IDRAIYCQ
Source	E.coli
Target Names	Os04g0533000
Protein Names	Recommended name: DEAD-box ATP-dependent RNA helicase 6 EC= 3.6.4.13
Expression Region	1-498
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.