



# Recombinant *Oryza sativa* subsp. japonica DEAD-box ATP-dependent RNA helicase 6 (Os04g0533000, LOC\_Os04g45040)

<b>Product Code</b>	CSB-MP773520OFG
<b>Abbreviation</b>	Os04g0533000, LOC_Os04g45040
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7XMK8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDPRARYPPG IGNGRGGNPN YYGRGPPPSQ HQQHQHQHQQ PPHPHHHQYV QRQPQPQQT HNSQHQQWLR RNQIAAEAAG ASEQKAPPVA DGIDSSSQDW KAQLKLPPQD TRYRTEDVTA TKGNEFEDYF LKRELLMGIY EKGFERPSPI QEESIPIALT GSDILARAKN GTGKTAAFCI PALEKIDQDK NAIQVVILVP TRELALQTSQ VCKELGKHLK IQVMVTTGGT SLKDDIVRLY QPVHLLVGTP GRILDLTCKG VCVLKNCSML VMDEADKLLS PEFQPSIQEL IRYLPSNRQI LMFSATFPVT VKEFKDKYLP KPYVINLMDE LTLKGITQFY AFVEERQKVH CLNTLFSKLQ INQSIIFCNS VNRVELLAKK ITELGYSCFY IHAKMLQDHR NRVFHDFRNG ACRNLVCTDL FTRGIDIQAV NVVINDFPK SAETYLHRVG RSGRFGHLGL AVNLITYEDR FNLYRIEQEL GTEIKPIPPQ IDRAIYCQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Os04g0533000
<b>Protein Names</b>	Recommended name: DEAD-box ATP-dependent RNA helicase 6 EC= 3.6.4.13
<b>Expression Region</b>	1-498
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .