



Recombinant *Oryza sativa* subsp. japonica Squamosa promoter-binding-like protein 8 (SPL8)

| | |
|--------------------------|--|
| Product Code | CSB-BP773522OFG |
| Abbreviation | SPL8 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q7XPY1 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. japonica (Rice) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MMNVPSAAAA SSCDDFGYNA TPPPPSLLP IMDQDGGGGS IQRDHHQHNN HQQLGYNLEP SSLALLPPSN AAAAAAHAT IAHASPHDLL QFYPTSHYLA AAGGAGGGGN PYSHFTAAAA AGSTFQSYQQ QPPQDAPEYY FPTLVSSAEE NMAFAATQL GLNLGYRTYF PPRGGYTYGH HPPRCQAEGC KADLSSAKRY HRRHKVCEHH SKAPVVVTAG GLHQRFCCQC SRFHLLDEFD DAKKSCRKRL ADHNRRRRKS KPSDGEHSGE KRRAQANKSA ATKDKAGSSS KNAGIGDGF TQLLGGAHMS KDQDQAMD LG EVVKEAVDPK GKASMQQQQQ QAHHGIHQQS HQQHGFPFPS SSGSCLFPQS QGAVSSTDTS NIAQVQEPSL AFHQQHHQHS NILQLGQAMF DLDFDH |
| Source | Baculovirus |
| Target Names | SPL8 |
| Protein Names | Recommended name: Squamosa promoter-binding-like protein 8 Alternative name(s): OsLG1 Protein LIGULELESS 1 |
| Expression Region | 1-416 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.