



# Recombinant *Oryza sativa* subsp. japonica Tubby-like F-box protein 7 (TULP7)

<b>Product Code</b>	CSB-BP773527OFG
<b>Abbreviation</b>	TULP7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7XSV4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSFRSIVRDF RDSFGTLSKR SFEVKISGFS GRHRGKSIGP SSELDDTPVV AQQSKWAGLP PELLRDVMKR LEEDDSNWPS RKDVVACASV CTTWRDMCKD IVRNPEFCGK LTFPVSLKQP GPRDGVICF IKRDKSKLTY HLYLCLSSAV LDETGKFLLS AKRSRRTTHT DYIISMDSKN ISRSSSGYIG KLRSNFLGTK FIIYDTQPPY NARTLCSQER TSTRFSSRKV SPKVPTGCYP IVQVNYELNV LGTRGPRRMQ CAMHSIPASA VEPGGIVPGQ PKELLPRLFE ESFRSMATSF SKYSITDHST DFSSSRFSEF GGGALQGQEQ EQDGDVDNKE RPLVLRNKAP RWHEQLQCWC LNFRGRVTV SVKNFQLIAA APQPSSGAAS EPSQAGQAAQ QQTQPSQPSS SSSSSSSNHD TVILQFGKVA KDMFTMDYRY PLSAFQAF AI CLTSFDTKLA CE
<b>Source</b>	Baculovirus
<b>Target Names</b>	TULP7
<b>Protein Names</b>	Recommended name: Tubby-like F-box protein 7 Short name= OsTLP7 Alternative name(s): Tubby-like F-box protein 3 Short name= OsTLP3
<b>Expression Region</b>	1-462
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.