



Recombinant *Oryza sativa* subsp. *japonica* Tubby-like F-box protein 7 (TULP7)

Product Code	CSB-YP773527OFG
Abbreviation	TULP7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XSV4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSFRSIVRDF RDSFGTLSKR SFEVKISGFS GRHRGKSIGP SSELDDTPVV AQQSKWAGLP PELLRDVMKR LEEDDSNWPS RKDVVACASV CTTWRDMCKD IVRNPEFCGK LTFPVSLKQP GPRDGVICF IKRDKSKLTY HLYLCLSSAV LDETGKFLLS AKRSRRTTHT DYIISMDSKN ISRSSSGYIG KLRSNFLGTK FIIYDTQPPY NARTLCSQER TSTRFSSRKV SPKVPTGCYP IVQVNYELNV LGTRGPRRMQ CAMHSIPASA VEPGGIVPGQ PKELLPRLFE ESFRSMATSF SKYSITDHST DFSSSRFSEF GGGALQGQEQ EQDGDVDNKE RPLVLRNKAP RWHEQLQCWC LNFRGRVTV SVKNFQLIAA APQPSSGAAS EPSQAGQAAQ QQTQPSQPSS SSSSSSSNHD TVILQFGKVA KDMFTMDYRY PLSAFQAFI CLTSFDTKLA CE
Source	Yeast
Target Names	TULP7
Protein Names	Recommended name: Tubby-like F-box protein 7 Short name= OsTLP7 Alternative name(s): Tubby-like F-box protein 3 Short name= OsTLP3
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.