



Recombinant Human Complement C1q-like protein 2 (C1QL2)

Product Code	CSB-EP773597HU-B
Abbreviation	C1QL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z5L3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	HYEMMGTCR MICDPYTAAP GGEPPGAKAQ PPGPSTAAL VMQDLSANPP PPFIQGPKGD PGRPGKPGPR GPPGEPGPPG PRGPPGEKGD SGRPGLPGLQ LTAGTASGVG VVGGGAGVGG DSEGEVTSAL SATFSGPKIA FYVGLKSPHE GYEVLKFDDV VTNLGNHYDP TTGKFSCQVR GIYFFTYHIL MRGGDGTSMW ADLCKNGQVR ASAIAQDADQ NYDYASNSVV LHLDSGDEVY VKLDGGKAHG GNNNKYSTFS GFLLYPD
Source	E.coli
Target Names	C1QL2
Protein Names	Recommended name: Complement C1q-like protein 2
Expression Region	22-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.