



Recombinant Human UPF0461 protein C5orf24 (C5orf24)

Product Code	CSB-EP773601HU-B
Abbreviation	C5orf24
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z6I8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMHPVASSNP AFCGPGKPCS LNEDAMRAAD QFDIYSSQQS KYSHTVNHKP MVCQRQDPLN ETHLQTTSGR SIEIKDELKK KKNLNRSGKR GRPSGTTKSA GYRTSTGRPL GTTKAAGFKT SPGRPLGTTK AAGYKVSPGR PPGSIKALSR LADLGYGCGT AAFPYPMMHG RAVHGVEETS SEVKPPNE
Source	E.coli
Target Names	C5orf24
Protein Names	Recommended name: UPF0461 protein C5orf24
Expression Region	1-188
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.