



# Recombinant Danio rerio Ankyrin repeat domain-containing protein 13C (ankrd13c)

<b>Product Code</b>	CSB-EP773623DIL
<b>Abbreviation</b>	ankrd13c
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7ZUV0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTGEKIRSVR KERKSGDLL EPDEEPAATG PAKHRGSKIF SGGNHRISRS SSSPGDPDGA YPVHECVFRG DVRRLSSLIR TQNIQKDVH GNTPLHLAVM MGHKECAHLL LAHNAPVKVK NAQGWSPLE AISYGDRQMI TALLRKLKQQ SRESVEDKRP RLLKALKELG DFYLELHWDF QSWVPLLSRI LPSDACKIYK QGINIRLDTT LIDFTDMKCQ RGDLSFIFCG DAPPESEFVV LDNEQKVYQR IHHEESEMET EEEVDILMSS DVYSATLSTK SITFSRAQTG WLFREDKTER VGNFLADFYM VNGLVLESRK RREHLSEEDI LRNKAIMESF SKGGSLIEQN FEPMRRQSLT PPSNTISWE EYITAETGKA PHLGRELKCK ESKKNFKATV AMSPDFPLGI ELLNVLEVI APFKHFNKLK EKVQMKLPPG FPKLDIPVF PTITATVTFQ EFRYDEFDES IFTIPSDYKG DPSRFPDL
<b>Source</b>	E.coli
<b>Target Names</b>	ankrd13c
<b>Protein Names</b>	Recommended name: Ankyrin repeat domain-containing protein 13C
<b>Expression Region</b>	1-488
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.