



# Recombinant Danio rerio S-methyl-5'-thioadenosine phosphorylase (mtap)

<b>Product Code</b>	CSB-YP773625DIL
<b>Abbreviation</b>	mtap
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7ZV22
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTSNSHVKIG IIGGSLDDP DILEGRTERY VVTPY GKPSD ALILGKIKNV DCVLLARHGR QHTIMPTNVN YQANIWALKE EGCTHLLVTT ACGSLREDIQ PGDIVLIDQF IDR TTKRVQT FYDGGQPTSP G VCHIPMAEP FCSKTREVL EVAQQLGVKC HTRGTMVTIE GPRFSSRAES LMFRQWGADV INMTTVPEVV LAKEAGLCYA SIAMATDYDC WKEHEEAVCV DNVLKT MKEN ANKASSILLT AIPQICQMDW DSTINAHKSM SQSSVMLPKH
<b>Source</b>	Yeast
<b>Target Names</b>	mtap
<b>Protein Names</b>	Recommended name: S-methyl-5'-thioadenosine phosphorylase EC= 2.4.2.28 Alternative name(s): 5'-methylthioadenosine phosphorylase Short name= MTA phosphorylase Short name= MTAP Short name= MTAPase
<b>Expression Region</b>	1-280
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.