



Recombinant Danio rerio Meteorin-like protein (metrnl)

Product Code	CSB-BP773626DIL
Abbreviation	metrnl
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZV46
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	QYSSDQCSW RGSGLTHEGH TRGVEQVYLR CAQGFLEWLY PTGAIIVNLR PNTLSPAASL LSVCIKPSKE SSGTHIYLDR LGKLRLLLSE GDQAEGKVHC FNIQDGALFI EAVPQRDISR KITAFQYELV NHRPGADPQS LSAPCQPCTD AEVLLAVCTS DFVARGRILG VSEEDEQTSV TVSLSHLYRQ KTQV FVSGGG RAKRWTGFVK MSRQCGVKPG DGEFLFTGTV RFGEAWLSCA PRYKDFLRVY QDARQQGTNP CHLETD
Source	Baculovirus
Target Names	metrnl
Protein Names	Recommended name: Meteorin-like protein
Expression Region	22-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.