



Recombinant Danio rerio Methionine-R-sulfoxide reductase B1-A (msrb1)

Product Code	CSB-EP773662DIL-B
Abbreviation	msrb1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q802G6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSFCSFSGGE IYKDFHESGM YVCAQCGYEL FSSRSKYEHS SPWPAFTETI HEDSVSKQEE RWGAYKVR CG KCGNGLGHEF VNDGPKHGLS RFUIFSSSLK FIPKVKNEQQ
Source	E.coli
Target Names	msrb1
Protein Names	Recommended name: Methionine-R-sulfoxide reductase B1-A Short name= MsrB Short name= MsrB1-A EC= 1.8.4.- Alternative name(s): Selenoprotein X-A Short name= SePR Short name= SelX-A
Expression Region	1-110
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.