



Recombinant Danio rerio UPF0415 protein C7orf25 homolog (zgc:55781)

Product Code	CSB-YP773665DIL
Abbreviation	zgc:55781
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q803H0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MATYNTLQDR IRVAKELLER VDKICSRQGR EVEGRAKLCG KLRAELKFLQ KVEAGKVVIK ESHLQSTNLT HLKAIVESAE SLEKVVSVLH VFAYESPDGQ KQTLVVDVVA NGGHTWVKS I GRKAEALHNI WQGRGQYGDK SVIRQAEDFL EASQQQPVQY SNPHIIFAFY NGVSSPMADK LKEMGISVRG DIVAVNTIEG GEEDEEDQE GDHEDLVEEE EDGEDDND DD SDDTDLMHTR VDRDTIVASL AFPTEVKVDV CNRVNLDITT LITYVSSL SH GNCHFTTFKEV VLTEQAAQER QEKVLPRL EE FMKGKELFAC QSAVEDFRVI LDTLGGPGEK SRAEELLARL KVVPDQP SER TQRLVMSSKV NRRSLMIFGT GDTLRAITMT ANSGFVRAAA NQGVRFSVFI HQPRALTEGK EWRATPI
Source	Yeast
Target Names	zgc:55781
Protein Names	Recommended name: UPF0415 protein C7orf25 homolog
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.